



LightStim ProPanel delivers a therapeutic, non-invasive treatment that provides clinically proven results to treat either anti-aging, acne, or pain. Relax, listen to music or meditate while enjoying this warm, soothing and restorative treatment.

Benefits Associated with LightStim ProPanel

- + Minimizes the appearance of pores, fine lines, and wrinkles
- + Helps energize and restore balance for optimal skin function
- + Promotes a more even skin tone
- + Calms and reduces redness
- + Produces a radiant glow
- + Works synergistically with other treatments and ingredients for better results
- + Increases blood circulation and oxygenation
- + Helps to accelerate recovery
- + Reduces and relieves pain
- + Destroys acne causing bacteria

ProPanel Anti-Aging Light

FDA Cleared to treat wrinkles on the entire face

ProPanel Acne Light

FDA Cleared to treat mild to moderate acne

ProPanel Pain Light

FDA Cleared to relieve arthritic pain, a variety of minor aches and pains, relax muscles, and increase local blood circulation.



Ingredient Guidelines

LightStim devices may be used on clean, dry skin or in conjunction with topical products. LightStim's unique technology works in harmony with high performance ingredients to enhance results.

- + Use LightStim PhotoSerum, LightStim SheetMasque or the appropriate products of your choice as described below.
- + When in doubt as to whether a topical product should be applied prior to or after LightStim treatment, ask the product manufacturer if it's safe to use a warm LED device on the skin after application of their product and/or simply apply product after LightStim treatment.
- + Ingredient categories that may be used prior to LightStim treatment (if non-irritating and safe to use with warmth) include collagen building peptides, cell stimulators, anti-inflammatories, antioxidants, and hydrators.
- + Ingredient categories to avoid prior to LightStim but that may be used after treatment include anything that may pose a risk of irritation or dryness, products that should not be exposed to warmth, Retin-A, benzoyl peroxide, harsh acids, potentially reactive forms of vitamin C and A, products containing metals (like colloidal silver and copper peptides), opaque masks and/or highly occlusive ingredients.



LightStim Precautions

- + Extreme sensitivity to light – migraines brought on by bright light
- + Use of photosensitizing drugs
- + Melasma or hyperpigmentation exacerbated by mild warmth
- + Do not apply products that should not be exposed to mild warmth prior to LightStim treatment
- + Accutane should be discontinued for 6 months in advance of treatment
- + Anyone who is pregnant must seek approval from a medical professional
- + Anyone with suspicious lesions or skin cancers should consult a medical professional
- + Anyone who is being treated for any type of serious health condition should seek approval from a medical professional before undergoing LED therapy

Instructions for Use

1. Plug into a wall outlet
2. Press the ON button
3. Set the desired number of minutes using the 5 MIN and 1 MIN increment buttons
 - Press 5 MIN 4 x for a standard 20-minute treatment
4. Press the SAVE button if you want to repeat this time increment in a future treatment
5. Press the ON button to turn the LEDs on and start the timer
6. The Light will automatically turn off at the end of the treatment time accompanied by a soft beep tone
7. To recall your saved time, press the RECALL button at any time





ProPanel Treatment Protocol

1. Cleanse
2. Apply treatment serums or transparent/near transparent mask if desired (refer to Ingredient Guidelines in this document)
3. Position LightStim provided goggles over the eyes and ensure the elastic band is in place snugly but comfortably around the head – if necessary moistened cotton rounds may be used under the goggles
4. Before moving the ProPanel into position over the face, turn the device on – this will provide an opportunity for the eyes to adjust to the brightness (see Instructions for Use above)
5. Position ProPanel over the face approximately 1- 2" away from the skin and adjust the 4 panels to conform to the face
6. When the treatment is complete, remove the device
7. Remove mask (if applied)
8. Apply moisturizer and sunscreen

Care & Cleaning

LEDs

You may clean the LEDs with alcohol once per day using a moderate amount lightly applied to a soft cloth. Gently wipe the LEDs. More than once per day may damage the LEDs over time. Do not use other cleaners on the LEDs except for alcohol as other cleaners/disinfectants will damage the LEDs and void your warranty.

Plastic Housing & Base

Use a gentle all-purpose cleaner, and if desired, your standard treatment room disinfectant cleaner. Spray the cleaner on a soft cloth and apply.

Do not spray alcohol or any other cleaners or disinfectants in or around the ProPanel. The spray can seep into the electronics, damage them, and void your warranty.

Do not soak the cloth in the alcohol/cleaner as the excess cleaner can seep into the electronics, damage them, and void your warranty.