

Saian Acne Protocol

Suggested Treatment time 60 minutes.

Suggested Price \$100.

1. Cleanse the neck and face with **Saian Active Cleanser**.
2. After thoroughly removing the cleanser, apply **Saian Salicylic Acid 2%** to the entire face with cotton rounds or glycolic applicator swabs, keeping away from the eyes and lips.
3. Remove the acid with cool water after 2 minutes.
4. Blend a teaspoon of **Saian Enzyme Mask** powder with equal part warm water to make a smooth paste, and apply to the face, neck and décolleté for 7-10 minutes with steam.
5. Remove the **Saian Enzyme Mask** thoroughly.
6. Spray the face with **Saian Herbal Toner** to help bring the pH back to normal.
7. Use a stiff fan brush to saturate the face with water, then cleanse the face with **Saian Professional Strength Ultrasonic Spatula** on the high setting (H) for 3-5 minutes. Make sure the skin is always wet while working with the device, and point the “blade” downwards at a 45 degree angle when in contact with skin.
8. Perform extractions.
9. Turn on your **Saian High Frequency machine**, and select the desired electrode for the treatment area. Perform the High Frequency treatment by whisking the wand gently over the face for 3-5 minutes. Gauze may be used between the skin and the electrode.
10. Apply a thin layer of **Saian Active Rx Mask** to the entire face for 10 minutes.
12. Remove the mask with warm water.
13. Apply 2 drops of **Saian Pure Hyaluronic Acid Serum** to the entire face.
14. After applying the serum, cover the client's eyes with protective goggles (optional for comfort) and position your **LightStim** device about an inch away from the client's face. Leave the mask under the light for 20 minutes for best results. The panel has a timer, and will shut off automatically.
15. Apply **Saian Acne Spot Treatment** using a q-tip to blemishes only. Do not shake the bottle prior to application.
16. Roll over the face with previously frozen **Saian Cool Duo Massagers** for 2-3 minutes to close the pores.
17. Spray the face with **Saian Sunscreen Mist**.

Suggested Home Care Products

>> **Saian Active Cleanser**

- Perfect for oily and blemished skin.

>> **Saian Hyaluronic Acid Serum**

- Hydrates the skin.

- Replenishes moisture without clogging the pores.

>> **Saian Acne Spot Treatment**

- Kills germs and bacteria and zaps blemishes overnight.

>> **Saian Sunscreen Mist**

- Water-based, oil-free, alcohol-free and fragrance-free.
- Excellent sun protection – SPF 30.
- Perfect for all skin types, especially formulated for acneic and sensitive skin.
- Perfect to re-apply after make-up.