

Saian FirmAge Non-Surgical Face-Lift

Suggested Treatment time 60-75 minutes (with/without facial massage)

Suggested Price \$250USD

1. Cleanse the neck and face with **Saian Replenishing Cleanser** or **Strawberry Cleanser**, depending on skin type.
2. Blend a teaspoon of **Saian Enzyme Mask** powder with equal part warm water to make a smooth paste, and apply to the face, neck and décolleté for 7-10 minutes with steam.
3. Remove the **Saian Enzyme Mask** thoroughly.
4. Spray the face with **Saian Herbal Toner** to help bring the pH back to normal.
5. Use a stiff fan brush to saturate the face with water, then cleanse the face with **Saian Professional Strength Ultrasonic Spatula** on the high setting (**H**) for 3-5 minutes. Make sure the skin is always wet while working with the device, and point the “blade” downwards at a 45 degree angle when in contact with skin.
6. Perform extractions and follow with **Saian High Frequency**, if needed.
7. Perform facial massage with **Essence by Saian Organic Beauty Oil**. (optional)
8. Remove the massage oil thoroughly.
9. Apply **Saian Hydrating Cellulosic Mask** according to instructions. Remove the mesh before applying the mask.
10. After applying the mask, cover the client’s eyes with protective goggles (optional for comfort) and position your **LightStim** anti-aging panel about an inch away from the client’s face. Leave the mask under the light for 15-20 minutes for best results. The panel has a timer, and will shut off automatically.
11. Remove the **Saian Hydrating Cellulosic Mask** thoroughly.
12. Apply a dime-sized amount of **Saian Neck and Chest Lift Crème** to the neck, jowls and décolleté.
15. Apply 1 pump of **Saian Active Renewal Serum** to the whole face focusing on wrinkle-prone areas.
16. Apply **Saian Active Renewal Crème** over the serum, focusing on wrinkle-prone areas.
17. Turn on the **Saian FirmAge Portable RF** device according to instructions. Set it to Low, Medium or High intensity. Use the device
26. Spray with Saian Sunscreen Mist SPF 30.

Suggested Retail Products

>> Saian Tri-Active Serum

- Protects and re-generates collagen, to firm and tone droopy and inelastic skin.
- Boosts synthesis of Collagen I, III, and IV.

- Uniforms fibril diameter and spacing.
- Inhibits enzymatic breakdown, avoiding excessive collagen damage.
- Has cryo-protective properties, which allows it to protect the skin's collagen even in the most severe cold environments.

>> **Saian Active Renewal Serum (best seller)**

- Improves wrinkle depth, wrinkle volume and roughness.
- Minimizes crow's feet and laugh lines.
- Prevents from further damage and wrinkle formation.

>> **Saian Active Renewal Crème (best seller)**

- Improves wrinkle depth, wrinkle volume and roughness.
- Minimizes crow's feet and laugh lines.
- Prevents from further damage and wrinkle formation.
- Hydrates the skin.

>> **SNC® by Saian Hyaluronic Acid Supplements** containing double strength hyaluronic acid with antioxidants for superior joint and spine lubrication, ease of mobility, skin and eye moisture, and overall increased elastin levels.