

Saian Ultimate Signature Facial

Suggested Treatment time 90 minutes

Suggested Price \$250USD

1. Cleanse the neck and face with **Saian Replenishing Cleanser** or **Strawberry Cleanser**, depending on skin type.
2. Blend a teaspoon of **Saian Enzyme Mask** powder with equal part warm water to make a smooth paste, and apply to the face, neck and décolleté for 7-10 minutes with steam.
3. Remove the **Saian Enzyme Mask** thoroughly.
4. Spray the face with **Saian Herbal Toner** to help bring the pH back to normal.
5. Use a stiff fan brush to saturate the face with water, then cleanse the face with **Saian Professional Strength Ultrasonic Spatula** on the high setting (**H**) for 3-5 minutes. Make sure the skin is always wet while working with the device, and point the “blade” downwards at a 45 degree angle when in contact with skin.
6. Perform extractions and follow with **Saian High Frequency**, if needed.
7. Perform facial massage with **Essence by Saian Organic Beauty Oil**. (optional)
8. Remove the massage oil thoroughly.
9. Apply **Saian Botanical Stem Cell Mask** to the skin in a thin layer avoiding the hairline and eyebrows/lashes. You will need approximately 2 pumps.
10. After applying the mask, cover the client’s eyes with protective goggles (optional for comfort) and position your **LightStim** anti-aging panel about an inch away from the client’s face. Leave the mask under the light for 15-20 minutes for best results. The panel has a timer, and will shut off automatically.
11. Remove the **Saian Botanical Stem Cell Mask** thoroughly.
12. After removing the mask, open the **Saian Hyaluronic Acid Ampoule** according to the directions on the box, making sure to wrap the top in a tissue before snapping it off. Pour the desired amount of serum into the palm of your hand, and apply to the face, neck and décolleté. Gently pat with open palms, especially around the wrinkle-prone areas.
13. Spray the face, neck and décolleté with **Saian Pure Collagen Elastin Spray** 4-6 times.
14. Apply a dime-sized amount of **Saian Neck and Chest Lift Crème** to the neck, jowls and décolleté.
15. Fold the **Saian Pure Collagen Sheet Mask** in half over the client’s face. Measure the distance from hairline to nostrils. Mark and cut the mask to expose the nostrils.
16. Apply the collagen mask immediately over the face.
17. Spray with **Saian Pure Collagen Elastin Spray** again over the mask 4-6 times just until it is moist.

18. Dip sponges into warm water, squeeze out excess water and press with moist sponges over the mask to mold it to the face, making sure there are no air bubbles.
19. Use the **Saian Skin Pro 2** portable galvanic on the “**nutrition**” mode rolling over the mask to increase absorption for 4 minutes. The machine has an automatic timer and will shut off after each 4 minute step.
20. Cover the mask with a warm towel, or turn on the steam over it for an additional 15 minutes.
21. Remove the collagen sheet by lifting it off in one piece.
22. Apply 2 drops of **Saian Tri-Active Collagen Serum** to the whole face focusing on wrinkle-prone areas.
23. Apply 2 drops of **Saian Pure Hyaluronic Acid Serum** to the whole face focusing on wrinkle-prone areas.
24. Apply **Saian Active Renewal Crème** (weightless, fragrance free) or **Saian Amino-Lift Leave-On Mask** (richer, antioxidant w/neroli oil) to maintain moisture.
25. Roll over the face with previously frozen **Saian Cool Duo Massagers** for 2-3 minutes to close the pores.
26. Spray with Saian Sunscreen Mist SPF 30.

Suggested Retail Products

>> Saian Tri-Active Serum

- Protects and re-generates collagen, to firm and tone droopy and inelastic skin.
- Boosts synthesis of Collagen I, III, and IV.
- Uniforms fibril diameter and spacing.
- Inhibits enzymatic breakdown, avoiding excessive collagen damage.
- Has cryo-protective properties, which allows it to protect the skin's collagen even in the most severe cold environments.

>> Saian Active Renewal Serum (best seller)

- Improves wrinkle depth, wrinkle volume and roughness.
- Minimizes crow's feet and laugh lines.
- Prevents from further damage and wrinkle formation.

>> Saian Active Renewal Crème (best seller)

- Improves wrinkle depth, wrinkle volume and roughness.
- Minimizes crow's feet and laugh lines.
- Prevents from further damage and wrinkle formation.
- Hydrates the skin.

>> SNC® by Saian Hyaluronic Acid Supplements containing double strength hyaluronic acid with antioxidants for superior joint and spine lubrication, ease of mobility, skin and eye moisture, and overall increased elastin levels.